Silk'n®

LED Light Therapy Mask

LumiLips FAC07NA

- EN USER MANUAL
- FR MODE D'EMPLOI
- ES MANUAL DEL USUARIO

Silk'n®

Device / Appareil / Dipositif


Table of contents

1. Warnings and precautions	2
2. Contraindications	7
3. Safety with Silk'n.....	8
3.1 Check your skin type.....	8
3.2 Acne severity.....	9
4. Get to know your device.....	10
4.1 Device description	10
4.2 Intended use	10
4.3 Package contents	10
4.4 Device parts	11
4.5 Device controls.....	11
4.6 Charging the battery	12
5. Treatment	13
5.1 Before first use.....	13
5.2 Treatment steps	14
6. Cleaning, maintenance & storage.....	15
7. Troubleshooting	16
8. Customer service.....	17
9. Replacement parts.....	17
10. Warranty	18
11. Disposal	18
12. Specifications	19
13. Environmental conditions.....	20
14. Signs and symbols	21
15. Labeling	23
16. Technical results.....	24
17. Legal notice	29

1. Warnings and precautions

This device is powerful and should be used with special attention to safety. Before first use, please read this user manual in full, and keep it for future use. Be sure to read all warnings, contraindications and safety precautions, and follow them strictly. We recommend you refamiliarize yourself with the instructions before each treatment.

The device is not intended for use by persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, without supervision or instruction concerning safe use of the equipment, to avoid the risk of fire or burns.

- Check if the voltage indicated on the device corresponds to the local mains voltage before use to avoid risk of electrocution or permanent damage to the device.
- If the wall socket used to charge the device has poor connections, the adapter becomes hot. Make sure you plug the device into a properly installed wall socket to avoid the risk of fire and burns.

- The adapter is not supplied with the device.
- If the charging cable is damaged, replace with a cable of the same standard (see section 9. Replacement parts).
- Always unplug the device after charging or if there is a power failure to avoid the risk of damage to the device.
- Do not use the device during charging.
- Do not touch the device, or charging cable with wet hands during charging.
- Never use the device near water or where water can be splashed onto the device, or if you think it has come into contact with water.
- The battery is non-replaceable. Do not modify or replace it.
- This device may cause radio frequency interference or may interfere with the operation of adjacent equipment. If interference occurs, appropriate mitigation measures need to be taken, such as repositioning or relocating the device or shielding the area.
- Do not look directly at the light source part of the device.

- Do not expose the device to harmful optical radiation.
- Before using the device you must carry out a sensitivity test (see section 5.1 Before first use).
- Only use the device as described in this user manual. Use of the device other than specified has not been tested and the risks are unknown.
- Do not use the device without removing creams, lotions or make-up from the face before treatment. Some ingredients in these products may cause light sensitivity or may reduce the effect of the treatment. If you use the device without removing these products, you may cause a skin reaction, or you may reduce the effect of the treatment.
- Do not use the device if the device, the charging cable or the adapter are damaged in any way.
- If the device malfunctions or operates abnormally, stop using it immediately and contact the manufacturer.
- If there is any change in the performance of the device, please stop using it immediately and contact the manufacturer in a timely manner.

- Turn off and remove the device immediately if it feels hot or your skin feels uncomfortably warm.
- Do not leave the device unattended when it is switched on to avoid the risk of fire or burns.
- Keep the device away from hot surfaces.
- Do not store flammable and explosive materials (sprayers, oil, gas, etc.) near the device.
- Do not use the device with heating accessories or near heat sources.
- Do not subject the device to heavy shocks to avoid risk of damage.
- Do not excessively bend or squeeze the charging cable.
- Do not expose the device to the outdoors for a long time to avoid contamination of the device by pests.
- Do not attempt to modify the device in any way. Modifying the device will void your warranty. Modifying the device could cause an electric shock when the device is used again. The device contains no serviceable parts.
- Do not connect this device with other devices not described in this user manual.

- Do not use accessories, removable parts or materials not described in this user manual.
- The device is not intended for use by children. Children should be supervised to ensure that they do not play with the equipment to avoid the risk of fire and burns.
- Do not place the device in the reach of pets or children to avoid misuse or damage to the device.
- Keep cables out of the reach of children. Cables of excessive length may cause strangulation.
- Avoid inhalation or ingestion of small parts.
- Persons with potential allergies to the accessible materials used in the device should not use the device.
- Only use accessories and cables supplied or approved by Silk'n. Use of unauthorised accessories or cables may increase electromagnetic emissions from the device.
- The device is a personal item and it is not recommended to share it with others.

2. Contraindications

Do not use the device to treat any other conditions apart from those listed in this user manual. The device has not been tested for any other conditions than those listed, and the risks are unknown.

Do not use if you:

- are pregnant, planning on becoming pregnant or breastfeeding. The device has not been tested on pregnant or breastfeeding women and therefore the risk to the foetus, neonate or pregnant women is unknown.
- have Fitzpatrick Skin Type I or Type V or Type VI (prone to hyper pigmentation).
- have a medical history of seizures triggered by light.
- have a history of photosensitive epilepsy.
- suffer from Lupus erythematosus, photosensitive eczema or Albinism. If you use the device to treat Lupus erythematosus, photosensitive eczema or Albinism you may cause a severe skin reaction.
- suffer from any photosensitive disorder (sensitivity to light). If you use the device and you suffer from a photosensitive disorder you may cause a severe skin reaction.
- are taking any medication that can cause photosensitivity. If you use the device and you are taking any medication that can cause photosensitivity you may cause a severe skin reaction. Photosensitivity is a common side effect of various medications. These can include certain antibiotics, chemotherapy drugs, and diuretics. If you are unsure about any medication you are taking, consult your healthcare provider.

- have severe acne.
- have had sunburn in the past 2 weeks.
- have open wounds (avoid treatment area until healed).

3. Safety with Silk'n

The device has the following safety features:

- **Auto shut-down:** The device automatically shuts down after 5 minutes of operation.
- **Low battery warning:** When the battery is low during operation, the buzzer will sound 5 times and the indicator light will flash red. The device will complete the current treatment cycle and then shut down. Charge the device before next use.

3.1 Check your skin type

Fitzpatrick skin types

The device is not suitable for all skin types. Use the table to compare your skin color to those in the table. This will help you determine if the device is suitable for you. Refer to the skin color chart at the back of this user manual.

Skin Type	Description	Can I use this device?
I	Very white skin, easily sunburned, never tanned.	X Do not use.
II	Fair skin, easy to burn, lightly tanned.	✓ Use with caution. Reduce the frequency of use.
III	Light olive skin, possible sunburn, gradually tanning.	✓ Use according to standard instructions.

Skin Type	Description	Can I use this device?
IV	Olive/light brown skin, easy to burn and tan.	✓ Use according to standard instructions.
V	Brown skin, rarely sunburned, easily tanned.	X Do not use.
VI	Dark brown/black skin, never sunburned.	X Do not use.

3.2 Acne severity

The device is not suitable for all acne conditions. Refer to the acne severity images at the back of this user manual to determine if the device is suitable for you.

Acne Severity	Description	Can I use this device?
Mild (I)	Some blackheads and whiteheads, a few pimples and pustules, no nodules.	✓ Yes
Moderate (II)	Many blackheads and whiteheads, some pimples and pustules, only a few small nodules.	✓ Yes
Severe (III)	Many blackheads, whiteheads, pimples and pustules, some nodules.	X No. This product is not recommended for severe acne as the benefits have not been established.

4. Get to know your device

4.1 Device description

The device is a lip treatment device which utilizes LED light therapy technology. The device combines 4 types of light therapy to treat the perioral region to reduce the appearance of fine lines and wrinkles and increase circulation.

The device includes LED lights, a mouthpiece, and a control button with indicator light. It has a rechargeable battery for cordless use.

It is an over-the-counter home use device for personal use. Do not use for any purpose other than as described in this user manual.

4.2 Intended use

The device is intended for the treatment of fine lines and wrinkles and to increase circulation within the perioral region.

4.3 Package contents

- LED lip mask
- USB Charging cable *
- User manual

* A Class II DC 5V, 0.5A mains adapter required to charge the device is not supplied.

4.4 Device parts

Refer to the diagram at the front of the manual to see the location of each feature.

- LED lip mask (#1)
- Control button with battery indicator light (#2)
- Charging port (#3)
- Mouthpiece (#4)
- LED lights (#5)
- USB charging cable (#6)

4.5 Device controls

Refer to the diagram at the front of the manual to see the location of each control.


Control	Action	Function / Mode
Control button with battery indicator light (#2) 	Long press (2 seconds).	The device starts and beeps once. The battery indicator light is white. The device operates continuously for 5 minutes and shuts down automatically.
	Long press (2 seconds) when the device is on.	The device beeps twice then shuts down.

Control	Action	Function / Mode
Battery / indicator light settings (#2)	White constant	Device on and ready for use.
	Green flashing	Battery charging
	Green constant	Battery charged
Battery / indicator light settings (#2)	Red flashing	Low battery
	During operation, a buzzer will sound 5 times and the indicator light will flash red.	The treatment will continue for the remaining 5 minutes and then shut down. You must charge the device before next use.
Charging port (#3) 		Port for the USB cable to charge the device.

4.6 Charging the battery

1. Fully charge the device before first use.
2. Connect the smaller end of the USB cable (#6) into the charging port (#3) of the device.
3. Connect the larger end of the USB cable (#6) into a compatible USB power source (e.g., a USB power adapter).
4. Connect the power adapter to the mains supply. The control button/battery indicator light (#2) flashes green.

- When the device is fully charged, the control button / battery indicator light (#2) is constant green.

-  You cannot use the device when charging the battery.
-  The battery icon does not change color.

5. Treatment

5.1 Before first use

Perform a sensitivity test before your first treatment. If you are sensitive to silicone, you should not use the device.


- Place the device over your inner arm.
 - Press and hold the control button (#2) for 2 seconds to start the device. Allow the device to run for the full 5-minute treatment program.
 - After the treatment test, wait at least 6 hours. If the skin on your inner arm appears red and/or itchy, the device is unsuitable for you. Return the device to the retailer within 14 days of purchase.
-  If your skin feels hot or uncomfortable, stop using the device immediately.

5.2 Treatment steps

⚠ Do not use the device without removing creams, lotions or makeup from the face before treatment. Some ingredients in these products may cause light sensitivity or may reduce the effect of the treatment.

⚠ Ensure that all parts of the device are properly assembled before use.

1. Insert the mouthpiece (**#4**) into the opening of the device (**#1**).
2. Ensure the teeth ridges and circular stem face upwards.

3. Push the mouthpiece forward until the circular stem passes through the opening and sits securely inside the device.
4. Clean and thoroughly dry your face and hands.
5. Insert the mouthpiece (**#4**) comfortably in your mouth.
6. Start the device. Long press the control button (**#2**) for 2 seconds. The device beeps once, and the indicator light (**#2**) turns on.
7. Relax during the 5-minute treatment. You may lie down if you wish.
8. The device operates continuously for 5 minutes, then shuts down automatically.

9. To shut down the device before the treatment ends, long press the control button (**#2**) for 2 seconds. The device beeps twice.
10. Remove the device and apply your preferred skincare serums or creams to moisturize the skin.

 If your skin feels hot or uncomfortable, remove the device from the skin. The device shuts down automatically after 5 minutes.

6. Cleaning, maintenance & storage

Cleaning

Clean the device after each treatment.

1. Use a soft, lint-free cloth dampened with water to wipe the device and the mouthpiece (**#4**).
 2. For stubborn residue, wipe with a soft cloth lightly dampened with up to 70% isopropyl alcohol.
 3. Ensure the device is completely dry before storing and next use.
-  Never use harsh detergents, acetone, or other abrasive liquids, because they will damage your device.
-  Never immerse your device or any of its parts in water.
-  Turn off the device before cleaning.
-  Do not clean any parts in the dishwasher.

Maintenance

- Do not perform maintenance while the device is in use or charging.
- The device contains no serviceable parts.
- Do not attempt to open, modify, or repair the device yourself.
- The battery is non-replaceable. Do not attempt to replace it.

Storage

- Store the device at room temperature in dry conditions.
- Avoid storing in dusty conditions, in direct sunlight, or in damp and humid conditions.

7. Troubleshooting

Do not attempt to open or repair your device. You may expose yourself to dangerous electrical components and may cause serious harm to yourself. It will also void your warranty.

Problem	Cause	Action
Cannot turn on the device.	The battery is low.	Charge the device.
LED lights do not light up.	The device is faulty.	Contact customer service.
Indicator light flashes red during use.	The battery voltage is too low.	Charge the device. The device will complete the current treatment cycle and then shut down.

Problem	Cause	Action
Device does not respond to button presses.	The device is faulty.	Contact customer service.

If a problem reoccurs and your device is within the warranty period, return it to your provider.

 Do not attempt to open or repair the device. If problems persist, contact customer service.

8. Customer service

For more information about Silk'n products, visit your regional Silk'n website: www.silkn.com or www.silkn.ca.

If the device is damaged, defective, in need of repair or if you need our assistance, please contact your nearest Silk'n service center. This manual is also available as a PDF download from the websites.

Service number: 1-800-250-5053 (toll-free)

Email: customercarena@silkn.com

9. Replacement parts

You may replace the 5V 0.5A charging cable with a cable of the same standard.

10. Warranty

EN

This product is covered by a 1-year warranty. The extent of the warranty on this product is limited to technical defects caused by faulty production processes. In the event you want to claim warranty, please contact our customer service for instructions. They may be able to solve your problem without needing to return the product to the store or our service center.

Our customer service is always pleased to help you!

- i** The manufacturer accepts no liability for damage or injuries caused by improper or incorrect use.

11. Disposal

Do not dispose of this device with household waste. This device contains electrical components and a lithium-ion battery. Dispose of the device in accordance with applicable local, state, or provincial regulations for electrical and electronic equipment.

Contact your local waste management authority for information on recycling and disposal options. Make sure the device is in an uncharged state before disposal.

Dispose of the packaging separated into single type materials. Dispose of cardboard and carton as waste paper and film via the recyclable material collection service.

12. Specifications

Model no.	LumiLips FAC07NA
Treatment area	Perioral region (area around the lips)
Technology	LED light therapy
LED light (wavelength; optical output power)	605nm + 630nm + 660nm + 880nm 40mW/cm ² ±20%
Operation & safety	Auto shutdown after 5 minutes
Dimensions	112 x 85 x 58 mm
System weight	86.5 gm
Risk group (optical radiation)	Exempt group
Operating conditions	Temperature: 5°C to 35°C
	Relative humidity: ≤85% RH
	Atmospheric pressure: 80 - 106 kPa
	Altitude: ≤2000 meters
Transport & storage between uses and storage condition	Temperature: -10°C to 50°C
	Relative humidity: ≤85% RH
	Atmospheric pressure: 80 - 106 kPa
Power input	DC 5V, 0.5A
Protection class	IP25
The device was tested and complies with these standards	IEC 60601-1 (Safety) IEC 60601-1-2 (EMC) IEC 60601-1-11 (Environment) IEC60601-2-83 (Safety of home light therapy equipment)

Expected service life of the equipment	1 year
Rechargeable Li-ion polymer battery	Charging time: Approximately 3 hours
	Nominal voltage: 3.7V
	Rated capacity: 500mAh, 1.85Wh
	Battery cycle life: ≥ 300 cycles

13. Environmental conditions

The device can withstand the environmental conditions for testing and transport listed in the table.

Temperature	$60 \pm 2^{\circ}\text{C}$ and $85\% \pm 2\%$ RH for 48 hours, then room temperature for 2 hours
	$-20 \pm 2^{\circ}\text{C}$ for 48 hours, then room temperature for 2 hours
	$40^{\circ}\text{C} \pm 3^{\circ}\text{C}$ and $90\% \pm 2\%$ RH for 48 hours, then room temperature for 2 hours
	40°C for 2 hours, with device functions tested in this environment
	0°C for 2 hours, with device functions tested in this environment

- i** Store the device for at least 48 hours before transfer to an operational environment and before use.

14. Signs and symbols

	A statement or event that indicates a possible or impending danger, problem when using the device.
	Important information about using or maintaining the device.
	For indoor use only.
IP25	IP classification. The device is protected against the ingress of objects larger than 12.5mm diameter and against water spray.
	Keep dry.
	Keep away from fire.
	Place upwards.
	Direct current.
	Do not dispose with household waste. Dispose in accordance with your country's laws and help protect the environment.
	Lithium-ion inside.

	The built-in rechargeable battery contains substances that may pollute the environment. Discard in accordance with local regulations.
	Interference may occur in the vicinity of the equipment.
	Degree of protection against electric shock: type BF applied part.
	Date of manufacture.
	Manufacturer.
	Product serial number.
	Product reference number.
	Medical device.
	Read the user manual.

15. Labeling

EN

- i** This device is not suitable for use in the presence of flammable anaesthetic mixture with air or with Oxygen or Nitrous Oxide.

* The image is for illustration purposes only.

16. Technical results

1. The LED Light Therapy Mask needs to take special precautions for electromagnetic compatibility (EMC), and the purchaser or user must install and use it according to the information about EMC Compliance provided in the instruction manual.
2. Portable and mobile radio frequency communication equipment may affect the LED Light Therapy Mask. Please keep away from such equipment when using it.
3. In addition to the connecting wires and charging cables sold by our company as spare parts of internal components, the use of accessories and cables outside the regulations may result in an increase in the emission of LED Light Therapy Mask or a decrease in immunity.

The following cable types must be used to ensure compliance with and adherence to interference radiation and immunity standards:

Cable	Length(m)
power cord	/
charging cable	1.0

4. This product should not be used close to or stacked with other equipment. If it must be used close or stacked, it should be observed to verify that it can operate normally under the configuration used.

5. This product should not touch the pins of the connectors marked with the electrostatic discharge warning symbol, and should not be connected to these connectors unless electrostatic discharge precautions are used.

Provisions on preventive measures against electrostatic discharge:

Human body or objects carry different electrostatic voltages due to charge transfer; since electrostatic discharge is completed in the order of ns or μ s, the peak current can reach tens of amperes, and the instantaneous power is very huge. The resulting electrostatic discharge electromagnetic pulse energy is enough to damage the sensitive components in the electronic components; because the rise time of the current waveform is very short, that is, the rate of change of the current (di/dt) is very large, it can induce a high potential of hundreds of volts or even thousands of volts, thereby a strong electric field is generated to break down the sensitive element. In order to prevent damage to the instrument, the following measures should be taken:

- (1) Ensure the environmental humidity.
- (2) Laying anti-static floor or carpet.
- (3) The operator should wear an anti-static hand strap on the wrist. This hand strap should have good grounding performance. This measure is the most effective.

6. The basic performance of the instrument:

Mode	Light	Radiation intensity $\pm 20\%$
1	Wavelength: 605nm+630nm+660nm+880nm	40mW/cm ²

7. The device has been tested and homologated in accordance with EN 60601-1-2 for EMC. This does not guarantee in any way that this device will not be affected by electromagnetic interference. Avoid using the device in high electromagnetic environment.

8. Do not use in MRI, diathermy, electrocautery, wireless power transfer, 5G, security systems and other environments with electromagnetic transmitters.

Technical Description Concerning Electromagnetic Emission		
Table 1: Declaration - electromagnetic emissions		
Guidelines and manufacturer's declaration-electromagnetic emissions		
The LED Light Therapy Mask is intended for use in the electromagnetic environment specified below. The customer or the user of the mask should ensure that it is used in such an environment.		
Launch test	Compliance	Electromagnetic environment-guide
RF emission CISRP11	Level of exemption	The LED Light Therapy Mask uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISRP11	Class B	The LED Light Therapy Mask is suitable for used in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emission IEC61000-3-2	Class A	
Voltage fluctuation/ flicker emission IEC61000-3-3	Complies	

Technical Description Concerning Electromagnetic Immunity			
Table 2: Guidance & Declaration - electromagnetic immunity			
Guidelines and manufacturer's declaration-electromagnetic immunity			
The LED Light Therapy Mask is intended for use in the electromagnetic environment specified below. The customer or the user of the mask should assure that It is used in such an environment.			
Immunity test	IEC60601 test level	Coincidence level	Electromagnetic environment-guide
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2, ±4, ±8, ±15kV air	±8kV contact ±2, ±4, ±8, ±15kV air	The floor should be wood, concrete, or ceramic tiles. If the floor is covered with a synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for the power lines ±1kV To input/ output lines	±2kV for the power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	±0.5, ±1kV line to line ±0.5, ±1, ±2kV line to earth	±0.5, ±1kV line to line ±0.5, ±1, ±2kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95% dip in UT.) for 0.5 cycle <5 % UT (>95% dip in UT.) for 1 cycle 70% UT (30% dip in UT) for 25 cycles <5% UT (>95 % dip in UT) for 250 cycles	<5 % UT (>95% dip in UT.) for 0.5 cycle <5 % UT (>95% dip in UT.) for 1 cycle 70% UT (30% dip in UT) for 25 cycles <5% UT (>95 % dip in UT) for 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the LED Light Therapy Mask requires continued operation during power mains interruptions, it is recommended that the LED Light Therapy Mask be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: UT is the AC mains voltage prior to application of the test level.			

Guidance & Declaration - electromagnetic immunity concerning Conducted RF & Radiated RF

Table 3 Guidance & Declaration - Electromagnetic immunity

The LED Light Therapy Mask is intended for use in the electromagnetic environment specified below. The customer or the user of the LED Light Therapy Mask should assure that it is used in such an environment.

Test	Test Level	Level	Environment Guide
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V 6V	Portable and mobile RF communications equipment should be used no closer to any part of the LED Light Therapy Mask, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \times P^{1/2}$ $d = 2 \times P^{1/2}$ $d = 1.2 \times P^{1/2}$ 80 MHz to 800 MHz $d = 2.3 \times P^{1/2}$ 800 MHz to 2.7 GHz P — the maximum rated output power of the transmitter provided by the transmitter manufacturer, in watts (W); d — the recommended separation distance, in meters (m). Field strengths from fixed RF transmitters are determined by an electromagnetic site survey "a", and in each frequency range "d" should be lower than the compliance level. Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted RF IEC 61000-4-6	6 Vrms ISM frequency band	3V/m	
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz		
	GHz		

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strengths from fixed transmitters, such as APP stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the LED Light Therapy Mask is used exceeds the applicable RF compliance level above, the LED Light Therapy Mask should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the LED Light Therapy Mask.

Table 4: Recommended separation distances between portable and mobile RF communications equipment and the LED Light Therapy Mask.

Recommended isolation distance between portable and mobile RF communication equipment and the LED Light Therapy Mask.

The LED Light Therapy Mask is intended for use in electromagnetic environment in which radiated RF disturbance is controlled. The customer or the user of the LED Light Therapy Mask can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the LED Light Therapy Mask as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150KHZ-80MHz $d=1.2p^{1/2}$	80MHz-800MHz $d=1.2p^{1/2}$	800MHz-2,7GHz $d=2.3p^{1/2}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For the maximum rated output power of the transmitter not listed in the table above, recommended separation distance "d" in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where "P" is the maximum rated output power of the transmitter provided by the transmitter manufacturer, in watts (W).

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

17. Legal notice

Copyright 2026 © Silkn Beauty Ltd. All rights reserved.

EN

Silkn Beauty Ltd. reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. Information furnished by Silkn Beauty Ltd. is believed to be accurate and reliable at the time of publication. However, Silkn Beauty Ltd. assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Silkn Beauty Ltd. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Silkn Beauty Ltd. Data is subject to change without notification. Silkn Beauty Ltd. has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property rights except as expressly provided in any written agreement from Silkn Beauty Ltd. Specifications are subject to change without notice.

Silk'n and the Silk'n logo are registered trademarks of Silkn Beauty Ltd., Galileistraat 17, 3029 AL Rotterdam, The Netherlands, www.silkn.com, www.silkn.ca, info@silkn.com.

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Copyright 2026 © Silkn Beauty Ltd. Tous droits réservés.
Copyright 2026 © Silkn Beauty Ltd. Todos los derechos reservados.
Silk'n Beauty Ltd.